
AC Adapter HHP-BFH01 Instruction Manual

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EN

Intended Use

This device is an AC Adapter intended for use with applicable Omron devices. The Omron devices which are compatible have this model mentioned on the rating label. In case compatibility is not clear please consult an authorized Omron service provider.

Please read this instruction manual and applicable Omron device manual thoroughly before using the unit. Please keep for future reference.

Important safety information

Warning: Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

- Do not use the AC adapter if the device or the power cord is damaged. Turn off the power and unplug the power cord immediately.
- Plug the AC adapter into the appropriate voltage outlet. Do not use in a multi-outlet plug.
- Never plug in or unplug the power cord from the electric outlet with wet hands.
- Do not disassemble or attempt to repair the AC adapter.

Caution: Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the equipment or other property.

- Fully insert the power plug into the outlet.
 - When disconnecting the power plug from the outlet, be sure to safely pull from the power plug.
 - Do not pull from the power cord.
 - When handling the power cord, take care not to do the following:
 - Do not damage. Do not break it.
 - Do not tamper with it. Do not forcibly bend or pull.
 - Do not twist. Do not bundle during use.
 - Do not pinch. Do not place under heavy objects.
 - Wipe any dust off of the power plug.
 - Unplug the monitor when not in use.
 - Disconnect the power plug before cleaning.
 - AC adapter cable may cause accidental strangulation in infants. Use this AC adapter only for compatible main Omron devices. Use of unsupported adapters may damage and/or may be hazardous to the device.
- While using this device, make sure that no mobile phone or any other electrical devices that emit electromagnetic fields is within 30cm. This may result in degradation of performance of the device.

Using the AC Adapter

Note: Make sure to use an easily accessible power socket in which to connect and disconnect the AC adapter.

- Insert the AC adapter plug into the AC adapter jack of the device.
- Plug the AC adapter into an electrical outlet.

To disconnect the AC adapter, unplug the AC adapter from the electrical outlet first, and then remove the AC adapter plug from the device.

Maintenance and Storage

To protect your device from damage, please follow the directions below:

- Store the device in a clean, safe location.
- Do not use any abrasive or volatile cleaners.
- Do not wash the device or immerse them in water.
- Do not use gasoline, thinners or similar solvents to clean the device.
- Disconnect the device from the electrical outlet and use a soft dry cloth to clean.
- Changes or modification not approved by the manufacturer will void the user warranty.

Do not disassemble or attempt to repair the device or components.

Specifications

Product Description	AC Adapter
Model	HHP-BFH01
Electrical Rating	AC100-240 V, 50-60 Hz, 0.12-0.065 A
Output voltage	6V DC / 700 mA
Durable period (Service life)	5 years
Operation conditions:	
Ambient temperature	10 to 40 °C
Relative humidity	15 to 90 % (non-condensing)
Ambient Pressure	700-1060 hPa

Transportation and storage conditions:

Ambient temperature	-20 to 60 °C
Relative humidity	10 to 90 % (non-condensing)
Weight	Approx. 67.2 g
Dimensions (box)	54 (w) × 72 (h) × 96 (l) mm
IP Classification*	IP22
Protection against electric shock	Class II ME equipment
Contents	AC Adapter, Instruction Manual

Notes:

- These specifications are subject to change without notice.
- IP classification is degrees of protection provided by enclosures in accordance with IEC 60529. The AC adapter is protected against solid foreign objects of 12.5 mm diameter and greater such as a finger. The AC adapter is protected against oblique falling water drops which may cause issues during a normal operation.
- This device can be used for continuous operation.
- Please report to the manufacturer and the competent authority of the Member State in which you are established about any serious incident that has occurred in relation to this device.

Warranty

Thank you for buying an OMRON product. This product is constructed of high quality materials and great care has been taken in its manufacturing. It is designed to give you every satisfaction, provided that it is properly operated and maintained as described in the instruction manual. This product is guaranteed by OMRON for a period of 1 year after the date of purchase. The proper construction, workmanship and materials of this product is guaranteed by OMRON. During this period of guarantee OMRON will, without charge for labour or parts, repair or replace the defect product or any defective parts. Should guarantee service be required, please contact your authorized OMRON distributor in your country. www.omron-healthcare.com Repair or replacement under the guarantee does not give rise to any extension or renewal of the guarantee period. The guarantee will be granted only if the complete product is returned together with the original invoice / cash ticket issued.

Important information regarding Electro Magnetic Compatibility (EMC)

HHP-BFH01 manufactured by OMRON HEALTHCARE Co., Ltd. conforms to EN60601-1-2:2015 Electro Magnetic Compatibility (EMC) standard. Further documentation in accordance with this EMC standard is available at OMRON HEALTHCARE EUROPE at the address mentioned in this instruction manual or at www.omron-healthcare.com. Refer to the EMC information for HHP-BFH01 on the website.

Correct Disposal of This Product (Waste Electrical & Electronic Equipment)

This marking shown on the product or its literature, indicates that it should not be disposed of, with other household wastes at the end of its working life. To prevent possible harm to the environment or human health from uncontrolled waste disposal, please separate this product from other types of wastes and recycle it responsibly to promote the sustainable reuse of material resources. Household users should contact either the retailer where they purchased this product, or their local government office, for details of where and how they can return this item for environmentally safe recycling. Business users should contact their supplier and check the terms and conditions of the purchase contract. This product should not be mixed with other commercial waste for disposal.

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