

OMRON

Comfort Cuff

Comfort Cuff

-
- **Instruction Manual**
 - **Mode d'emploi**
 - **Gebrauchsanweisung**
 - **Manuale di istruzioni**
 - **Manual de instrucciones**
 - **Gebruiksaanwijzing**
 - **РУКОВОДСТВО ПО ЭКСПЛУАТАЦИИ**
-

EN

FR

DE

IT

ES

NL

RU

Issue Date:

Date de publication :

Ausgabedatum:

Data di pubblicazione:

Fecha de publicación: Uitgiftedatum:

2019-10-15

Дата выпуска:

OMRON


Comfort Cuff

Instruction Manual for OMRON's Comfort Cuff

Thank you for purchasing the OMRON Comfort Cuff. Arm circumference for this cuff is 22 to 42 centimeters. OMRON's Comfort cuff can be connected to various OMRON Blood Pressure Monitors. Please refer to the original instruction manual of your Blood Pressure Monitor.

Please follow this instruction manual thoroughly for your safety. Please keep for future reference. For specific information about your own bloodpressure, CONSULT YOUR PHYSICIAN

Packaging Contents



- A. Arm cuff (arm circumference 22-42 cm)
- B. Air tube
- C. Air plug
- D. Air plug


Instructions for use:

Note: Connect correct air plug to the tube according to your device.


Applying the Arm Cuff

You can wrap the cuff either on your right or left arm.

- Remove tight-fitting clothing from your upper arm.
- Do not place the cuff over thick clothes and do not roll up your sleeve if it is too tight.




1. Insert the air plug into the air jack (on the left side of the device).




The cuff must be fully deflated when it is inserted into the air jack.

2. Sit in a chair with your feet flat on the floor and place your arm on a table so that the cuff will be at the same level as your heart.




3. Hold the grip on the cuff securely with your hand.



4. Turn the palm of your hand upward.



5. Turn the palm of your hand upward.



Apply the cuff to your upper arm so that the blue strip is centred on the middle of your inner arm and points down the inside of the arm. The air tube should run down the inside of your forearm and be in line with your middle finger. The bottom of the cuff should be approximately 1 to 2 cm above the elbow.

6. When the cuff is positioned correctly, close the fabric fastener **FIRMLY**.




Notes:

- Make certain the cuff fits snugly around your arm.
- The cuff should make good contact with your skin. You should be able to fit your index finger between the cuff and your arm easily, so you can pull the cuff off and on.
- Make sure that there are no kinks in the air tubing.

Taking measurements on the right arm

Please note the following points when applying the cuff to your right arm. Apply the cuff so that the air tube is at the side of your elbow.



Notes:

- Be careful not to rest your arm on the air tube, or otherwise restrict the flow of air to the cuff.
- Apply the arm cuff so that no part of the cuff is positioned over the elbow joint. The cuff should be 1 to 2 cm above the elbow.
- If used 2 times per day, the cuff will last for 3 years.

Maintenance

Stains on the cuff can be carefully removed using a moist cloth and soapy water.

Do not wash the cuff or make it wet!
















Do not fold the cuff or air tube together too tightly and do not crease them.







EN

Specifications:

Product description	OMRON Blood Pressure Monitor Cuff
Model	Comfort Cuff
Arm circumference	22 to 42 cm
Cuff dimensions	Approx. 151 mm × 582 mm
Cuff pressure range	Pressure: 0 to 299 mmHg
Accuracy	Pressure: ± 3 mmHg
Classification	Class I
Operating conditions	+10 to +40°C (50 to 104°F) 15 to 90% RH (non-condensing) 700 to 1060 hPa
Storage/Transport conditions	-20 to +60°C (-4 to 140°F) 10 to 95% RH (non-condensing) 700 to 1060 hPa
Package contents	Cuff, instruction manual, air plug
Durable period (Service life), Warranty	1 year

Symbols description

	CE Marking		Medical device
	LOT number		Humidity limitation
	Temperature limitation		Identifier of cuffs compatible for the device
	Atmospheric pressure limitation	INDEX 	Range pointer and brachial artery alignment position
	Cuff positioning indicator for the left arm	LATEX FREE	Not made with natural rubber latex
	Artery mark		Reference catalogue number
	Range indicator of arm circumferences to help selection of the correct cuff size		Intelli Wrap Cuff
	Easy Cuff		Comfort Cuff

	Technology and Design, JAPAN		Date of manufacture
	Technology by OMRON Healthcare in Japan		Technology and Quality, JAPAN
	Arm circumference		Need for the user to consult the instruction manual
Product production date is integrated in the Lot number, which is placed on the product and/or sales package: the first 4 digits mean year of production, the next 2 digits mean month of production and the last 2 digits mean day of production.			



Please report to the manufacturer and the competent authority of the Member State in which you are established about any serious incident that has occurred in relation to this device.

Manufacturer 	OMRON HEALTHCARE Co., Ltd. 53, Kunotsubo, Terado-cho, Muko, KYOTO, 617-0002 JAPAN
EU-representative 	OMRON HEALTHCARE EUROPE B.V. Scorpius 33, 2132 LR Hoofddorp, THE NETHERLANDS www.omron-healthcare.com
Importer in EU	
Production facility	OMRON DALIAN Co., Ltd. No. 3, Song Jiang Road, Economic and Technical Development Zone, Dalian 116600, China
Subsidiaries	OMRON HEALTHCARE UK LTD. Opal Drive, Fox Milne, Milton Keynes, MK15 0DG, UK www.omron-healthcare.com
	OMRON MEDIZINTECHNIK HANDELSGESELLSCHAFT mbH Konrad-Zuse-Ring 28, 68163 Mannheim, GERMANY www.omron-healthcare.com
	OMRON SANTÉ FRANCE SAS 3, Parvis de la Gare, 94130 Nogent-sur-Marne, FRANCE www.omron-healthcare.com

EN