OMRON

Automatic Upper Arm Blood Pressure Monitor

M2+ (HEM-7188-LE)


M2 Essential (HEM-7188-E)

English

Instruction Manual 1: Safety and Other Information

For symbols information, refer to “Symbols Description” section.

Intelli sense

Introduction

Thank you for purchasing the OMRON Automatic Upper Arm Blood Pressure Monitor. This blood pressure monitor uses the oscillometric method of blood pressure measurement. This means this monitor detects your blood movement through your brachial artery and converts the movements into a digital reading.

Safety Instructions
This instruction manual provides you with important information about the OMRON Automatic Upper Arm Blood Pressure Monitor. To ensure the safe and proper use of this monitor, READ and UNDERSTAND all of these instructions. **If you do not understand these instrons or have any questions, contact your OMRON retail outlet or distributor before attempting to use this monitor. For specific information about your own blood pressure, consult with your physician.**

Intended Use

Intended Purpose
This device is a digital monitor intended for use in measuring blood pressure and pulse rate in adult patient population. The device detects the appearance of irregular heartbeats during measurement and indicates this via a symbol with the measurement result.

Intended Patients

Adult patient population
AC Adapter (optional accessory) Handling and Usage
Fully insert the AC adapter into the outlet.

- When unplugging the AC adapter from the outlet, be sure to safely pull on the AC adapter. DO NOT pull on the AC adapter cable.
- DO NOT touch it./ DO NOT break it./ DO NOT tamper with it./ DO NOT pinch it./ DO NOT forcibly bend or pull it.
- Read and follow the “Correct Disposal of this Product” in section 5 when disposing of the device and any used accessories or optional parts.

AC Adapter (optional accessory) Handling and Usage
Fully insert the AC adapter into the outlet.

- When unplugging the AC adapter from the outlet, be sure to safely pull on the AC adapter. DO NOT pull on the AC adapter cable.
- DO NOT touch it./ DO NOT break it./ DO NOT tamper with it./ DO NOT pinch it./ DO NOT forcibly bend or pull it.
- Read and follow the “Correct Disposal of this Product” in section 5 when disposing of the device and any used accessories or optional parts.

Clinical Benefit
Patient’s blood pressure can be measured non-invasively and simply in the home environment.
Type of Use
This monitor is intended to be multiple patient multiple use.

Limitation
Patient’s arm circumference must be 17 – 42 cm.
Indication
This device is used by healthy individuals, patients with hypertension, patients with health-conscious individuals, in a general household or office for measuring blood pressure and pulse rate.

Receiving and Inspection
Remove this monitor and other components from the packaging and inspect for damage. If this monitor or any other components is damaged, DO NOT USE and consult with your OMRON retail outlet or distributor.

Important Safety Information

Read the Important Safety Information in this instruction manual before using this monitor. Follow this instruction manual thoroughly for your safety.
For specific information about your own blood pressure, consult with your physician.

Contraindications
DO NOT use this monitor on an injured arm or an arm under medical treatment.
DO NOT use the arm cuff on your arm while on an intravenous drip or blood transfusion.
DO NOT use this monitor on infants, toddlers, children or persons who are unable to sense themselves.

Side Effects

- Taking measurements more often than necessary may cause bruising due to blood flow interference.
- Inflating to a higher pressure than necessary may result in bruising of the arm where the cuff is applied. NOTE: Refer to “If your systolic pressure is more than 210 mmHg” in the end of instruction manual 2 for additional information.
- Stop using this monitor and consult with your physician if you experience skin irritation or discomfort.

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

- DO NOT adjust medication (including altering your use of any drug or treatment) based on readings from this blood pressure monitor. Take medication as prescribed by your physician. ONLY a physician is qualified to diagnose and treat high blood pressure and heart related conditions.
- NEVER diagnose or treat yourself based on your readings. ALWAYS consult with your physician.
- If you are experiencing any symptoms or concerns, consult with your physician.
- DO NOT postpone/stop regular checkups or physician visits based on the results that you are getting from this monitor.
- DO NOT use this monitor in areas containing high frequency (HF) surgical equipment, magnetic resonance imaging (MRI) equipment, computerized tomography (CT) scanners. This may result in incorrect operation of the monitor and/or cause an inaccurate reading.
- DO NOT use this monitor in oxygen rich environments or near flammable gas.
- Consult with your physician before using this monitor if you have common arrhythmias such as atrial or ventricular premature beats or atrial fibrillation; arterial sclerosis; poor perfusion; diabetes; pregnancy; pre-eclampsia or renal disease. NOTE that any of these conditions in addition to patient motion, trembling, or shivering may affect the measurement reading.
- To help avoid strangulation, keep the air tube and AC adapter cable away from infants, toddlers and children.
- This product contains small parts that may cause a choking hazard if swallowed by infants, toddlers and children.



AC Adapter (optional accessory) Handling and Usage
DO NOT use the AC adapter if this monitor or the AC adapter cable is damaged. If this monitor or the cable is damaged, turn off the power and unplug the AC adapter immediately.
Plug the AC adapter into the appropriate voltage outlet. DO NOT use the multi-outlet plug.
NEVER plug in or unplug the AC adapter from the electric outlet with wet hands.
DO NOT disassemble or attempt to repair the AC adapter.

Battery Handling and Usage

- Keep the batteries out of the reach of infants, toddlers and children.

Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient, or cause damage to the equipment or other property.

- Consult with your physician before using this monitor on an arm where intravascular access or therapy, or an antero-venous (A-V) fistula is present because of temporary interference to blood flow and could result in injury.
- Consult with your physician before using this monitor if you have any mastectomy or lymph node clearance.
- Consult with your physician before using this monitor if you have severe blood flow problems or blood disorders as cuff inflation can cause bruising.
- ONLY inflate the arm cuff when it is applied on your upper arm.
- Remove the arm cuff if it does not start deflating during a measurement.
- DO NOT use this monitor for any purpose other than measuring blood pressure.
- During measurement, make sure that no mobile device or any other electrical device that emits electromagnetic fields is within 30 cm of this monitor. This may result in incorrect operation of the monitor and/or cause an inaccurate reading.
- DO NOT disassemble or attempt to repair this monitor or other components. This may cause an inaccurate reading.
- DO NOT use in a location where there is moisture or a risk of water splashing this monitor. This may damage this monitor.
- DO NOT use this monitor in a moving vehicle such as in a car or on an aircraft.
- DO NOT drop or subject this monitor to strong shocks or vibrations.
- DO NOT use this monitor in places with high or low humidity or high or low temperatures. Refer to section 4.
- During measurement, observe the arm to ensure that the monitor is not causing prolonged impairment to blood circulation.
- DO NOT use this monitor with other medical electrical equipment simultaneously. This may result in incorrect operation of the devices and/or cause an inaccurate reading.
- Avoid bathing, drinking alcohol or caffeine, smoking, exercising and eating for at least 30 minutes before taking a measurement.
- Rest for at least 5 minutes before taking a measurement.
- Remove tight-fitting or cloth clothing from your arm while taking a measurement.
- Remain still and DO NOT talk while taking a measurement.

Display/Problem	Possible Cause	Solution
 appears		Apply the arm cuff correctly. Then take another measurement. Refer to section 4 of instruction manual 2. Remain still and sit correctly during a measurement. If the irregular heartbeat symbol
 appears	The pulse rate is not detected correctly.	 continues to appear. It is recommended that you consult with your physician.
 does not flash during a measurement		 continues to appear. It is recommended that you consult with your physician.
 flashes	Batteries are low.	Replacing all 4 batteries with new ones is recommended. Refer to section 3 of instruction manual 2.
 and  appears, or the monitor is turned off unexpectedly during a measurement.	Batteries are depleted.	Immediately replace all 4 batteries with new ones. Refer to section 3 of instruction manual 2.
Nothing appears on the display of the monitor.	Battery polarities are not properly aligned.	Check the battery installation for proper placement. Refer to section 3 of instruction manual 2.
Readings appear too high or too low.	Blood pressure varies constantly. Many factors including stress, time of day, and/or how you apply the arm cuff, may affect your blood pressure. Review sections 2, 4 and 5 of instruction manual 2.	Press the [START/STOP] button to turn the monitor off, then press it again to take a measurement. If the problem continues, remove all batteries and wait for 30 seconds. Then re-install the batteries. If the problem still persists, contact your OMRON retail outlet or distributor.
Any other problem occurs.		Press the [START/STOP] button to turn the monitor off, then press it again to take a measurement. If the problem continues, remove all batteries and wait for 30 seconds. Then re-install the batteries. If the problem still persists, contact your OMRON retail outlet or distributor.

2 Limited Warranty

Thank you for buying an OMRON product. This product is constructed of high quality materials and great care has been taken in its manufacturing. It is designed to give you a high level of comfort, provided that it is properly operated and maintained as described in the instruction manual.
Household users should be warned by OMRON for a period of 5 years after the date of purchase. The proper construction, workmanship and materials of this product is warranted by OMRON. During this period of warranty OMRON will, without charge for labour or parts, repair or replace the defect product or any defective parts. The warranty does not cover any of the following:
A. Transport costs and/or maintenance.
B. Costs for repairs and/or defects resulting from repairs done by unauthorised persons.
C. The device check-ups and maintenance.
D. Failure or wear of optional parts or other attachments other than the main device itself, unless explicitly warranted above.
E. Costs arising due to non-acceptance of a claim (those which will be charged for).
F. Damages of any kind including personal accident accidentally or from misuse.
G. Calibration service is not included within the warranty.
H. Optional parts have one (1) year warranty from date of purchase.
Optional parts include, but are not limited to the following items:

- Cuff and cuff tube.

Should warranty service be required please apply to the dealer whom the product was purchased from or an authorised OMRON distributor. For the address refer to the product packaging / literature or to your specialised retailer. If you have difficulties in finding OMRON customer services, visit our website (www.omron-healthcare.com) for contact information.
Repair or replacement under the warranty does not give rise to any extension or renewal of the warranty period.
The warranty will be granted only if the complete product is returned together with the original invoice/ cash ticket issued to the consumer by the retailer.

Should warranty service be required please apply to the dealer whom the product was purchased from or an authorised OMRON distributor. For the address refer to the product packaging / literature or to your specialised retailer. If you have difficulties in finding OMRON customer services, visit our website (www.omron-healthcare.com) for contact information.

AC Adapter (optional accessory) Handling and Usage
Fully insert the AC adapter into the outlet.

- When unplugging the AC adapter from the outlet, be sure to safely pull on the AC adapter. DO NOT pull on the AC adapter cable.
- DO NOT touch it./ DO NOT break it./ DO NOT tamper with it./ DO NOT pinch it./ DO NOT forcibly bend or pull it.
- Read and follow the “Correct Disposal of this Product” in section 5 when disposing of the device and any used accessories or optional parts.

Battery Handling and Usage

- DO NOT insert batteries with their polarities incorrectly aligned.
- ONLY use 4 “AA” alkaline or manganese batteries with this monitor. DO NOT use other types of batteries. DO NOT use new and used batteries together. DO NOT use different brands of batteries together.
- Remove the batteries if this monitor will not be used for a long period of time.
- If battery fluid should get in your eyes, immediately rinse with plenty of clean water. Consult with your physician immediately.
- If battery fluid should get on your skin, wash your skin immediately with plenty of fluid, lukewarm water. If irritation, injury or pain persists, consult with your physician.
- DO NOT use batteries after their expiration date.
- Periodically check the batteries to ensure they are in good working condition.

General Notices

- To stop a measurement, press the [START/STOP] button while taking a measurement.
- ONLY use 4 “AA” alkaline or manganese batteries with this monitor. DO NOT use other types of batteries. DO NOT use new and used batteries together. DO NOT use different brands of batteries together.
- Periodically check the batteries to ensure they are in good working condition.

⚠️ Precaution
DO NOT disassemble or attempt to repair this monitor or other components. This may cause an inaccurate reading.

3.2 Storage

- Store your monitor and other components in a clean, safe location.
- Gently fold the air tube into the arm cuff. Do not bend or crease the air tube excessively.
- Do not store your monitor and other components:
 - If your monitor and other components are wet.
 - In locations exposed to extreme temperatures, humidity, direct sunlight, dust or corrosive vapors such as bleach.
 - In locations exposed to vibrations or shocks.

3.3 Wiping the Monitor

- Do not use any abrasive or volatile cleaner.
- Use a soft dry cloth or a soft cloth moistened with mild (neutral) detergent to wipe your monitor and arm cuff, and then wipe them with a dry cloth.
- Do not wash or immerse your monitor and arm cuff or other components in water.
- Do not use gasoline, thinners or similar solvents to wipe your monitor and arm cuff or other components.

3.4 Calibration and Service

- The accuracy of this blood pressure monitor has been carefully tested and is designed for a long service life.
- It is generally recommended to have the unit inspected every two years to ensure correct functioning and accuracy. Please consult your authorised OMRON dealer or the OMRON Customer Service at the address given on the packaging or attached literature.

4 Specifications	
Product Category	Electronic Sphygmomanometers
Product description	Automatic Upper Arm Blood Pressure Monitor
Model (Code)	M2+ (HEM-7188-LE) / M2 Essential (HEM-7188-E)
Display	LCD digital display
Cuff pressure range	0 to 299 mmHg
Blood pressure measurement range	SYS: 60 to 260 mmHg <p>DIA: 40 to 115 mmHg</p>
Pulse measurement range	40 to 180 beats / min.
Accuracy	Pressure: ±3 mmHg <p>Pulse: ±5 % of display reading</p>
Measurement method	Oscillometric method
Operation mode	Continuous operation
IP classification	Monitor: IP21 <p>Optional AC adapter: IP21 (HHP-AM01 / HHP-AM01) or IP22 (HHP-BF01)</p>
Rating	DC6 V 4 W
Power source	4 “AA” batteries 1.5 V or optional AC adapter (INPUT AC 100 - 240 V 50/60 Hz 0.12 - 0.065 A)
Battery life	Approx. 900 measurements (Using new alkaline batteries and included arm cuff). Depending on the type of battery and arm cuff)
Durable period (Service life)	Monitor: 5 years or the time when reaching 30000 times of use. / Cuff: 5 years or the time when reaching 10000 times of use. / Optional AC adapter: 5 years
Operating conditions	+10 to +40 °C / 15 to 90 % RH (non-condensing) / 800 to 1060 hPa
Storage / Transport conditions	-20 to +60 °C / 10 to 90 % RH (non-condensing)
Weight	Monitor: approx. 260 g (not including batteries) <p>Arm cuff for HEM-7188-E: approx. 170 g <p>Arm cuff for HEM-7188-E: approx. 110 g</p></p>
Dimensions	Monitor: approx. 98 mm (W) x 79 mm (H) x 138 mm (L) <p>AC adapter: approx. 145 mm x 594 mm (air tube: 750 mm) <p>Arm cuff for HEM-7188-E: approx. 146 mm x 466 mm (air tube: 610 mm)</p></p>
Cuff circumference applicable to the monitor	17 to 42 cm (included arm cuff) <p>HEM-7188-LE: 22 to 42 cm <p>HEM-7188-E: 22 to 32 cm)</p></p>

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
Intelli sense

3.1 Maintenance

- To protect your monitor from damage, follow the directions below. Changes or modifications not approved by the manufacturer will void the user warranty.

⚠️ Precaution
DO NOT disassemble or attempt to repair this monitor or other components. This may cause an inaccurate reading.

3.2 Storage

- Store your monitor and other components in a clean, safe location.
- Gently fold the air tube into the arm cuff. Do not bend or crease the air tube excessively.
- Do not store your monitor and other components:
 - If your monitor and other components are wet.
 - In locations exposed to extreme temperatures, humidity, direct sunlight, dust or corrosive vapors such as bleach.
 - In locations exposed to vibrations or shocks.

3.3 Wiping the Monitor

- Do not use any abrasive or volatile cleaner.
- Use a soft dry cloth or a soft cloth moistened with mild (neutral) detergent to wipe your monitor and arm cuff, and then wipe them with a dry cloth.
- Do not wash or immerse your monitor and arm cuff or other components in water.
- Do not use gasoline, thinners or similar solvents to wipe your monitor and arm cuff or other components.

3.4 Calibration and Service

- The accuracy of this blood pressure monitor has been carefully tested and is designed for a long service life.
- It is generally recommended to have the unit inspected every two years to ensure correct functioning and accuracy. Please consult your authorised OMRON dealer or the OMRON Customer Service at the address given on the packaging or attached literature.

4 Specifications

Product Category Electronic Sphygmomanometers
Product description Automatic Upper Arm Blood Pressure Monitor
Model (Code) M2+ (HEM-7188-LE) / M2 Essential (HEM-7188-E)
Display LCD digital display
Cuff pressure range 0 to 299 mmHg
Blood pressure measurement range SYS: 60 to 260 mmHg
DIA: 40 to 115 mmHg
Pulse measurement range 40 to 180 beats / min.
Accuracy Pressure: ±3 mmHg
Pulse: ±5 % of display reading
Measurement method Oscillometric method
Operation mode Continuous operation
IP classification Monitor: IP21
Optional AC adapter: IP21 (HHP-AM01 / HHP-AM01) or IP22 (HHP-BF01)
Rating DC6 V 4 W
Power source 4 “AA” batteries 1.5 V or optional AC adapter (INPUT AC 100 - 240 V 50/60 Hz 0.12 - 0.065 A)
Battery life Approx. 900 measurements (Using new alkaline batteries and included arm cuff). Depending on the type of battery and arm cuff)
Durable period (Service life) Monitor: 5 years or the time when reaching 30000 times of use. / Cuff: 5 years or the time when reaching 10000 times of use. / Optional AC adapter: 5 years
Operating conditions +10 to +40 °C / 15 to 90 % RH (non-condensing) / 800 to 1060 hPa
Storage / Transport conditions -20 to +60 °C / 10 to 90 % RH (non-condensing)
Weight Monitor: approx. 260 g (not including batteries)
Arm cuff for HEM-7188-E: approx. 170 g
Arm cuff for HEM-7188-E: approx. 110 g
Dimensions Monitor: approx. 98 mm (W) x 79 mm (H) x 138 mm (L)
AC adapter: approx. 145 mm x 594 mm (air tube: 750 mm)
Arm cuff for HEM-7188-E: approx. 146 mm x 466 mm (air tube: 610 mm)
Cuff circumference applicable to the monitor 17 to 42 cm (included arm cuff)
HEM-7188-LE: 22 to 42 cm
HEM-7188-E: 22 to 32 cm)

3.1 Maintenance

- To protect your monitor from damage, follow the directions below. Changes or modifications not approved by the manufacturer will void the user warranty.

⚠️ Precaution
DO NOT disassemble or attempt to repair this monitor or other components. This may cause an inaccurate reading.

3.2 Storage

- Store your monitor and other components in a clean, safe location.
- Gently fold the air tube into the arm cuff. Do not bend or crease the air tube excessively.
- Do not store your monitor and other components:
 - If your monitor and other components are wet.
 - In locations exposed to extreme temperatures, humidity, direct sunlight, dust or corrosive vapors such as bleach.
 - In locations exposed to vibrations or shocks.

3.3 Wiping the Monitor

- Do not use any abrasive or volatile cleaner.
- Use a soft dry cloth or a soft cloth moistened with mild (neutral) detergent to wipe your monitor and arm cuff, and then wipe them with a dry cloth.
- Do not wash or immerse your monitor and arm cuff or other components in water.
- Do not use gasoline, thinners or similar solvents to wipe your monitor and arm cuff or other components.

3.4 Calibration and Service

- The accuracy of this blood pressure monitor has been carefully tested and is designed for a long service life.
- It is generally recommended to have the unit inspected every two years to ensure correct functioning and accuracy. Please consult your authorised OMRON dealer or the OMRON Customer Service at the address given on the packaging or attached literature.

4 Specifications

Product Category Electronic Sphygmomanometers
Product description Automatic Upper Arm Blood Pressure Monitor
Model (Code) M2+ (HEM-7188-LE) / M2 Essential (HEM-7188-E)
Display LCD digital display
Cuff pressure range 0 to 299 mmHg
Blood pressure measurement range SYS: 60 to 260 mmHg
DIA: 40 to 115 mmHg
Pulse measurement range 40 to 180 beats / min.
Accuracy Pressure: ±3 mmHg
Pulse: ±5 % of display reading
Measurement method Oscillometric method
Operation mode Continuous operation
IP classification Monitor: IP21
Optional AC adapter: IP21 (HHP-AM01 / HHP-AM01) or IP22 (HHP-BF01)
Rating DC6 V 4 W
Power source 4 “AA” batteries 1.5 V or optional AC adapter (INPUT AC 100 - 240 V 50/60 Hz 0.12 - 0.065 A)
Battery life Approx. 900 measurements (Using new alkaline batteries and included arm cuff). Depending on the type of battery and arm cuff)
Durable period (Service life) Monitor: 5 years or the time when reaching 30000 times of use. / Cuff: 5 years or the time when reaching 10000 times of use. / Optional AC adapter: 5 years
Operating conditions +10 to +40 °C / 15 to 90 % RH (non-condensing) / 800 to 1060 hPa
Storage / Transport conditions -20 to +60 °C / 10 to 90 % RH (non-condensing)
Weight Monitor: approx. 260 g (not including batteries)
Arm cuff for HEM-7188-E: approx. 170 g
Arm cuff for HEM-7188-E: approx. 110 g
Dimensions Monitor: approx. 98 mm (W) x 79 mm (H) x 138 mm (L)
AC adapter: approx. 145 mm x 594 mm (air tube: 750 mm)
Arm cuff for HEM-7188-E: approx. 146 mm x 466 mm (air tube: 610 mm)
Cuff circumference applicable to the monitor 17 to 42 cm (included arm cuff)
HEM-7188-LE: 22 to 42 cm
HEM-7188-E: 22 to 32 cm)

3.1 Maintenance

- To protect your monitor from damage, follow the directions below. Changes or modifications not approved by the manufacturer will void the user warranty.

⚠️ Precaution
DO NOT disassemble or attempt to repair this monitor or other components. This may cause an inaccurate reading.

3.2 Storage

- Store your monitor and other components in a clean, safe location.
- Gently fold the air tube into the arm cuff. Do not bend or crease the air tube excessively.
- Do not store your monitor and other components:
 - If your monitor and other components are wet.
 - In locations exposed to extreme temperatures, humidity, direct sunlight, dust or corrosive vapors such as bleach.
 - In locations exposed to vibrations or shocks.

3.3 Wiping the Monitor

- Do not use any abrasive or volatile cleaner.
- Use a soft dry cloth or a soft cloth moistened with mild (neutral) detergent to wipe your monitor and arm cuff, and then wipe them with a dry cloth.
- Do not wash or immerse your monitor and arm cuff or other components in water.
- Do not use gasoline, thinners or similar solvents to wipe your monitor and arm cuff or other components.

3.4 Calibration and Service

- The accuracy of this blood pressure monitor has been carefully tested and is designed for a long service life.
- It is generally recommended to have the unit inspected every two years to ensure correct functioning and accuracy. Please consult your authorised OMRON dealer or the OMRON Customer Service at the address given on the packaging or attached literature.

4 Specifications

Product Category Electronic Sphygmomanometers
Product description Automatic Upper Arm Blood Pressure Monitor
Model (Code) M2+ (HEM-7188-LE) / M2 Essential (HEM-7188-E)
Display LCD digital display
Cuff pressure range 0 to 299 mmHg
Blood pressure measurement range SYS: 60 to 260 mmHg
DIA: 40 to 115 mmHg
Pulse measurement range 40 to 180 beats / min.
Accuracy Pressure: ±3 mmHg
Pulse: ±5 % of display reading
Measurement method Oscillometric method
Operation mode Continuous operation
IP classification Monitor: IP21
Optional AC adapter: IP21 (HHP-AM01 / HHP-AM01) or IP22 (HHP-BF01)
Rating DC6 V 4 W
Power source 4 “AA” batteries 1.5 V or optional AC adapter (INPUT AC 100 - 240 V 50/60 Hz 0.12 - 0.065 A)
Battery life Approx. 900 measurements (Using new alkaline batteries and included arm cuff). Depending on the type of battery and arm cuff)
Durable period (Service life) Monitor: 5 years or the time when reaching 30000 times of use. / Cuff: 5 years or the time when reaching 10000 times of use. / Optional AC adapter: 5 years
Operating conditions +10 to +40 °C / 15 to 90 % RH (non-condensing) / 800 to 1060 hPa
Storage / Transport conditions -20 to +60 °C / 10 to 90 % RH (non-condensing)
Weight Monitor: approx. 260 g (not including batteries)
Arm cuff for HEM-7188-E: approx. 170 g
Arm cuff for HEM-7188-E: approx. 110 g
Dimensions Monitor: approx. 98 mm (W) x 79 mm (H) x 138 mm (L)
AC adapter: approx. 145 mm x 594 mm (air tube: 750 mm)
Arm cuff for HEM-7188-E: approx. 146 mm x 466 mm (air tube: 610 mm)
Cuff circumference applicable to the monitor 17 to 42 cm (included arm cuff)
HEM-7188-LE: 22 to 42 cm
HEM-7188-E: 22 to 32 cm)

3.1 Maintenance

- To protect your monitor from damage, follow the directions below. Changes or modifications not approved by the manufacturer will void the user warranty.

⚠️ Precaution
DO NOT disassemble or attempt to repair this monitor or other components. This may cause an inaccurate reading.

3.2 Storage

- Store your monitor and other components in a clean, safe location.
- Gently fold the air tube into the arm cuff. Do not bend or crease the air tube excessively.
- Do not store your monitor and other components:
 - If your monitor and other components are wet.
 - In locations exposed to extreme temperatures, humidity, direct sunlight, dust or corrosive vapors such as bleach.
 - In locations exposed to vibrations or shocks.

3.3 Wiping the Monitor

- Do not use any abrasive or volatile cleaner.
- Use a soft dry cloth or a soft cloth moistened with mild (neutral) detergent to wipe your monitor and arm cuff, and then wipe them with a dry cloth.
- Do not wash or immerse your monitor and arm cuff or other components in water.
- Do not use gasoline, thinners or similar solvents to wipe your monitor and arm cuff or other components.

3.4 Calibration and Service

- The accuracy of this blood pressure monitor has been carefully tested and is designed for a long service life.
- It is generally recommended to have the unit inspected every two years to ensure correct functioning and accuracy. Please consult your authorised OMRON dealer or the OMRON Customer Service at the address given on the packaging or attached literature.

4 Specifications

Product Category Electronic Sphygmomanometers
Product description Automatic Upper Arm Blood Pressure Monitor
Model (Code) M2+ (HEM-7188-LE) / M2 Essential (HEM-7188-E)
Display LCD digital display
Cuff pressure range 0 to 299 mmHg
Blood pressure measurement range SYS: 60 to 260 mmHg
DIA: 40 to 115 mmHg
Pulse measurement range 40 to 180 beats / min.
Accuracy Pressure: ±3 mmHg
Pulse: ±5 % of display reading
Measurement method Oscillometric method
Operation mode Continuous operation
IP classification Monitor: IP21
Optional AC adapter: IP21 (HHP-AM01 / HHP-AM01) or IP22 (HHP-BF01)
Rating DC6 V 4 W
Power source 4 “AA” batteries 1.5 V or optional AC adapter (INPUT AC 100 - 240 V 50/60 Hz 0.12 - 0.065 A)
Battery life Approx. 900 measurements (Using new alkaline batteries and included arm cuff). Depending on the type of battery and arm cuff)
Durable period (Service life) Monitor: 5 years or the time when reaching 30000 times of use. / Cuff: 5 years or the time when reaching 10000 times of use. / Optional AC adapter: 5 years
Operating conditions +10 to +40 °C / 15 to 90 % RH (non-condensing) / 800 to 1060 hPa
Storage / Transport conditions -20 to +60 °C / 10 to 90 % RH (non-condensing)
Weight Monitor: approx. 260 g (not including batteries)
Arm cuff for HEM-7188-E: approx. 170 g
Arm cuff for HEM-7188-E: approx. 110 g
Dimensions Monitor: approx. 98 mm (W) x 79 mm (H) x 138 mm (L)
AC adapter: approx. 145 mm x 594 mm (air tube: 750 mm)
Arm cuff for HEM-7188-E: approx. 146 mm x 466 mm (air tube: 610 mm)
Cuff circumference applicable to the monitor 17 to 42 cm (included arm cuff)
HEM-7188-LE: 22 to 42 cm
HEM-7188-E: 22 to 32 cm)

3.1 Maintenance

- To protect your monitor from damage, follow the directions below. Changes or modifications not approved by the manufacturer will void the user warranty.

⚠️ Precaution
DO NOT disassemble or attempt to repair this monitor or other components. This may cause an inaccurate reading.

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




Instruction Manual 2: Operational Instructions

- ET Kasutusjuhend 2: kasutussuunistid
- LV Lietošanas instrukcija 2: ekspluatācijas instrukcijas
- LT Naudojimo instrukcija Nr. 2: naudojimo instrukcijos
- BG Ръководство за употреба 2: Инструкции за работа
- SR Uputstvo za upotrebu 2: uputstva za rukovanje
- RO Manual cu instrucțiuni 2: Instrucțiuni de utilizare

Read Instruction Manual 1 and 2 before use.

Tutvuge enne kasutamist kasutusjuhenditega 1 ja 2.
 Pirms lietošanas izlasiet lietošanas instrukciju 1. un 2.
 Prieš naudodami, perskaitykite naudojimo instrukcijas Nr. 1 ir 2.
 Преди употреба прочетете ръководството за употреба 1 и 2.
 Pre upotrebe pročitajte uputstvo za upotrebu 1 i 2.
 Cititi manualul de instrucțiuni 1 și 2 înainte de utilizare.

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- A. Tube side of the cuff should be 1 - 2 cm above the inside of your elbow.
- B. Make sure that the air tube is on the inside of your arm and wrap the cuff securely so it can no longer slip round.

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- Relax and sit comfortably. Remain still and do not talk.
- 1. Keep your back and arm supported.
- 2. Keep the arm cuff at the same level as your heart.
- 3. Keep your feet flat and your legs uncrossed.

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- 1. Press the [START/STOP] button.
- A: Flashes at every heartbeat.
- B: Appears while the cuff is deflating.
- 2. The reading is saved automatically.

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