

Comparison Table

Digital and Infrared Thermometers



Model	GentleTemp	Flex Temp Smart	Eco Temp Smart	Eco Temp Basic	i-Temp	i-Temp mini	O-Temp III
Model Number	MC-510-E2	MC-343F-E	MC-341-E	MC-246-E	MC-670-E	MC-271W-E	MC-204-E2
Use	Ear	Oral, Axillar, Rectal	Oral, Axillar, Rectal	Oral, Axillar, Rectal	Oral, Axillar	Oral, Axillar, Rectal	Oral, Axillar, Rectal
Beeper	●	●	●	●	●	●	●
Memory	Last Measurement	Last Measurement	Last Measurement	Last Measurement	Last Measurement	Last Measurement	Last Measurement
Exchangeable Battery	●	●	●	●	●	●	●
Celsius & Fahrenheit	●	●	●	●	Celsius	●	●
3 Digit Display (1 decimal place)	●	●	●	●	●	●	●
Water resistant		●	●	●		●	●
Flexible Tip		●					
Automatic Shut Off	●	●	●	●	●	●	●